

PlumePort®

Laparoscopic Smoke Filtration Device

Improves Visibility

Safely filters surgical smoke improving surgeon's visibility, while protecting staff and patients from airborne hazards.

Adjustment Flow Control

Lightweight design and integral control maintains peritoneal distention.

Moisture Management

With an integrated fluid trap and wicking materials, moisture is effectively managed to ensure the integrity of the filtration media.

ULPA Filtration


Contains more than two times ULPA filtration media than competitive devices providing enhanced protection and extending filter life.

Activated Carbon

Incorporates activated carbon known to adsorb gases.

Universal Connectivity

Easily connects to all standard luer lock fittings.



Product Ordering Information

Description	Catalog Number	Quantity
PlumePort	PP010*	10/Case

*Only available outside of the U.S.

