

PlumePen[®] Pro

Surgical Smoke Evacuation Pencil

Sleek, Slim Profile

Smaller than competitive surgical smoke evacuation pencils; similar in size to standard ESU pencils

Removable Smoke Capture Port

Allows for easy electrode and capture port replacement

Compatible with Standard Electrodes

Designed for surgeries where longer electrodes are desired

Full Blade Tip Reveal

Provides an unobstructed view of surgical site

360° Swivel

Allows for fluid hand movement

Textured Profile

Improves surgeon's grip



Additional capture ports available to accommodate 4 in (102 mm), 5 in (127 mm) and 6 in (152 mm) electrodes.




PlumePen® Pro

Surgical Smoke Evacuation Pencil

The new, improved PlumePen® Pro features an easy to remove x-ray detectable capture port that allows for effortless electrode replacement to accommodate surgeries where longer electrodes are desired. Pro is smaller than competitive surgical smoke evacuation pencils and similar in size to standard electrosurgical pencils.

Pro is a cost effective and suitable alternative for surgeries where a pencil adapter is less optimal.



Product Ordering Information

Description	Catalog Number	Quantity
PlumePen® Pro (Non-Stick ESU Blade) 10 ft (3 m) tubing	PLPRO4020	20/Case
PlumePen® Pro (Non-Stick ESU Blade) 15 ft (4.6 m) tubing	PLPRO4520	20/Case
Capture Ports (X-ray detectable) (Accommodates 4 in [102 mm] electrodes)	PROTIP420	20/Case
Capture Ports (X-ray detectable) (Accommodates 5 in [127 mm] electrodes)	PROTIP520	20/Case
Capture Ports (X-ray detectable) (Accommodates 6 in [152 mm] electrodes)	PROTIP620	20/Case

Available Capture Ports



PROTIP420 . PROTIP520 . PROTIP620
(each sold separately)

