


PlumePort® SEO

Laparoscopic Smoke Filtration Device



Controlled Management

Manages smoke throughout the surgical procedure.

Universal Connection

Standard luer lock fitting; connects to most surgical smoke evacuators.

Up to 50% More Flow

Provides greater flow than passive laparoscopic devices.

Adjustable Flow Rate

Allows flow adjustment via surgical smoke evacuator.

Excellent Value

An economical effective solution

Quiet and Easy to Use

Technical Information — Flow Rates (Nominal Values)

ViroVac® on Medium (8.3 lpm)

PlumeSafe™ TURBO (60%) (5.1 lpm)

PlumeSafe™ TURBO (9.2 lpm) (40% in TURBO mode)

VisiClear® (Lapmode 90% = 6 lpm)

Product Ordering Information

Description	Catalog Number	Quantity
PlumePort® SEO	PP012SEO	12/Case

