

# LaparoVue™ Visibility System



*All-in-One Solution*  
*for improved laparoscopic visibility*


# Laparoscopic Visualization

## Complexity Made Simple



### warms & white balances

- warms scope in 5 minutes to prevent fogging
- unit warms up to 5 hours
- allows for white balancing




### cleans

- provides crisp, clear view in 5 seconds
- minimizes procedural delays
- radiopaque VueTip™ Trocar Swabs remove debris and excess liquid from trocar cannula
- x-ray detectable Microfiber Cloth designed to capture residue and particulates at a microscopic level

### defogs

- pre-filled warming solution helps prevent intra-operative scope fogging



\* Photographs are representative of a fogged and clear / defogged lens but do not reflect actual product use.

# Distinctive By Design

## Unique Product Attributes

### *one unit, two ports*

- two ports: warming / white balancing and cleaning / defogging
- fits 3mm–12mm scopes without port size converter
- allows for simultaneous use of both ports

### *detachable scope cradle*

- safely and securely supports scope
- reduces tension on lens and potential for damage

### *built-in LED lights*

- provides enhanced port visibility

### *ergonomic design*


- stable base prevents tipping
- contoured, easy grip profile
- engineered with soft silicone material



Actual product may differ slightly in appearance to images shown.


# Change Your Vue . . .

for laparoscopic & robotic-assisted laparoscopic surgeries



replace all this . . .

with LaparoVue™ Visibility System  
an All-in-One Solution



### Product Ordering Information

Description	Catalog Number	Quantity
LaparoVue™ Visibility System (includes main unit/scope cradle, 2 trocar swabs, and 1 microfiber cloth)	LAPVUE10	10/Case
VueTip™ Trocar Swab Accessory Kit (includes 2 trocar swabs and 1 microfiber cloth)	VUETIP10	10/Case



VueTip™ Trocar Swab Accessory Kit  
(sold separately; part: VUETIP10)



**Trial Now — call 1.716.835.7000 or 1.800.343.2324**

To view additional **Laparoscopic Solutions** visit [www.buffalofilter.com](http://www.buffalofilter.com).