

# ViroVac<sup>®</sup>

## Surgical Smoke Evacuator



### **Compact Design**

*Allows for easy placement on shelf or counter and accommodates stackability.*

### **Patented 3-Port Filter Design**

*Versatile 3-port filter design features covered, hidden male ports and accommodates a variety of tubing sizes providing versatility in procedure and accessory setup.*

### **Variable Filter Life (up to 35 hours)**

*Delivers unsurpassed value and reduces per procedure costs.*

### **Filter Tracking Technology**

*Automatically tracks and records filter life through a patented technology.*

### **Simple LED Display**

*Two button design (motor speed and on/standby) makes this system extremely easy to operate.*

### **Quiet Operation**

*Addresses noise concerns.*



# ViroVac®

## Surgical Smoke Evacuator

ViroVac surgical smoke plume evacuator provides ideal flow for physician offices and outpatient surgery centers. It combines customizable (variable) filter life, a compact shelf design, 3-port filter with filter life tracking technology, while exceeding expectations for value and efficacy.



### Technical Information

**Noise Levels** (Nominal Values)  
@ maximum power < 55 dBA

**Dimensions** (H x W x D)  
6 in x 11 in x 15.5 in  
15.2 cm x 27.9 cm x 39.4 cm

**Flow Rates**  
25 cfm (708 lpm) using 7/8 in (22 mm) tubing

**Weight**  
With filter installed: 12.0 lbs (5.0 kg)  
Without filter installed: 9.5 lbs (4.0 kg)

**Warranty**  
One (1) Year Full Warranty


### Product Ordering Information

Description	Catalog Number	Quantity
ViroVac® (100/120 Volt Systems) Includes one filter and one footswitch	VV120	1/Each
ViroVac® (200/240 Volt Systems) Includes one filter and one footswitch	VV220	1/Each
ViroVac® (100/120 Volt Systems) Includes one footswitch	VV120 NF	1/Each
V (200/240 Volt Systems) Includes one footswitch	VV220 NF	1/Each
EZLink® Automatic activation device, can be used to remotely turn the evacuator on and off, saving filter life.	EZLINK01	1/Case
ViroSafe® Fluid Trap Designed to capture fluids to help maintain the integrity of the filtration system and assist with maintaining filter life.	VSFT10	10/Case
ViroSafe® ViroVac® Filter	VS353	1/Case
Four layers of proven filtration gross particulate, ULPA media, activated carbon, and post filter.	VS35302	2/Case
	VS35304	4/Case

### Accessories (each sold separately)



EZLink®



ViroSafe® Fluid Trap

