

# PenAdapt<sup>®</sup>

## Surgical Smoke Evacuation Adapter

### Universal Design

*Pliable silicone material allows for proper fit with most electrosurgical pencil brands.*

### Angled Back End Clasp


*Reduces tension placed on surgeon's wrist providing comfort.*

### Efficient Workflow

*Eliminates need for OR personnel to hold tubing.*

### Low Profile Front Nozzle

*Improves visibility of the surgical site.*



PA1025, PA1025CB,  
PA1020, PA1020CB



PA2010B, PA2010BCB

#### Product Ordering Information

Description	Catalog Number	Quantity
PenAdapt Ideal for wall suction applications. Includes 3/8 in (9.5 mm) x 1 ft (30.5 cm) tubing with an adapter to 3/8 in (9.5 mm) tube.	PA1025	25/Case
PenAdapt (Cut-Back) Ideal for wall suction applications. Cut-back adapter provides increased visibility. Includes 3/8 in (9.5 mm) x 1 ft (30.5 cm) tubing with an adapter to 3/8 in (9.5 mm) tube.	PA1025CB	25/Case
PenAdapt Ideal for surgical smoke evacuator applications. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube.	PA2010	10/Case
PenAdapt Ideal for surgical smoke evacuator applications. Cut-back adapter provides increased visibility. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube.	PA2010CB	10/Case
PenAdapt Ideal for attachment to reusable electrosurgical pencil applications. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube.	PA2010B	10/Case
PenAdapt Ideal for attachment to reusable electrosurgical pencil applications. Cut-back adapter provides increased visibility. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube.	PA2010BCB	10/Case

\*ESU pencils sold separately.



# SnapEvac™

## Surgical Smoke Evacuation Adapter

### Snap and Evac

*Just snap in handpiece of your choice*

### Efficient Workflow

*Eliminates need for OR personnel to hold tubing*

### Transparent Capture Port

*Improves visibility of the surgical site*

### Cost-Effective

*Economical option for surgical smoke evacuation*

### Lightweight Tubing



#### Product Ordering Information

Description	Catalog Number	Quantity
SnapEvac™ Designed using ABS plastic and a polycarbonate capture port with 10 ft (3 m) tubing; 7/8 in (22 mm) cuff.	SnapEvac20	20/Case

\*ESU pencil and handpiece both sold separately. SnapEvac may not fit all ESU pencils and handpieces.

