

PlumeSafe® Turbo

Surgical Smoke Evacuator



(Filter Sold Separately)

Turbo Mode

Significantly increases the airflow to quickly remove smoke from the surgical site.

Interactive LCD Display

Puts information at your fingertips, multiple languages are included, adding to the unit's versatility.

Adjustable Flow Control

Suction setting can easily be adjusted in 10% increments.

Filter Tracking Technology

Automatically tracks and records filter life through a patented technology.

Quiet Operation

Addresses noise concerns.

Product Ordering Information

Description	Catalog Number	Quantity
PlumeSafe® Turbo (100/120 Volt 50/60 Hz) Includes one footswitch; filter sold separately	PSWTURBO AS	1/Case
PlumeSafe® Turbo (220/240 Volt 50/60 Hz) Includes one footswitch; filter sold separately	PSWTURBO AS/220	1/Case
EZLink Automatic activation device, can be used to remotely turn the evacuator on and off, saving filter life.	EZLINK01	1/Case

Technical Information

Noise Levels (Nominal Values)	@maximum power < 58 dBa
Dimensions (H x W x D)	13.5 in (34 cm) x 8 in (20 cm) x 17 in (43 cm)
Flow Rates	68 cfm (1923 lpm) using 1-1/4 in (32 mm) tubing 44 cfm (1257 lpm) using 7/8 in (22 mm) tubing
Weight	With filter installed: 16.9 lbs (7.0 kg) Without filter installed: 14.0 lbs (6.0 kg)
Warranty	Three (3) Year Full Warranty

Accessories (each sold separately)



EZLink®



BUFFALO FILTER®
Filtration Group®